

Induction validation - Cell Culture

Introduction

Get started by giving your protocol a name and editing this introduction.

Materials

- >
- >

12 well plates
- >

Dox
- >

TMP
- >

Media with puro (2 ug/mL)
- >

15 mL falcon tubes
- >

1.5 mL eppendorf tubes

Procedure

Conditions

1. Biological triplicates of the following conditons

Total = 30

Table1

	A	B	C	D	E
1		Dox (nM)	TMP(uM)	total media volume needed	total cells needed
2	CRISPRi				
3	GFI1B w/drugs	600	5		
4	GFI1B w/o drugs	0	0		
5	CD46 w/drugs	600	5		
6	CD46 w/o drugs	0	0		
7	NTC w/o drugs	0	0		
8	CRISPRa				
9	GFI1B w/drugs	600	0		
10	GFI1B w/o drugs	0	0		
11	CD46 w/drugs	600	0		
12	CD46 w/o drugs	0	0		
13	NTC w/o drugs	0	0		

Day 0

2. Warm media (RPMI with puro)

3. Count cells

CRISPRi GF11B =

CRISPRi CD46 =

CRISPRi NTC =

CRISPRa GF11B =

CRISPRa CD46 =

CRISPRa NTC =

4. Pellet 320,000 cells per condition

CRISPRi GF11B t =

CRISPRi GF11B ut =

CRISPRi CD46 t =

CRISPRi CD46 ut =

CRISPRi NTC ut =

CRISPRa GF11B t =

CRISPRa GF11B ut =

CRISPRa CD46 t =

CRISPRa CD46 ut =

CRISPRa NTC ut=

5. Make working dilutions of Dox and TMP

Add 1 uL of Dox stock (50 mg/mL) to 999 uL media/puro

Add 10 uL of TMP stock (100 mM) to 990 uL media/puro

6. Make media + drug mixes

Add 15 mL of media/puro with 92.4 uL of WD Dox and 75 uL WD TMP

Add 15 mL of media/puro with 92.4 uL of WD Dox

7. Resuspend with 6.4 mL of respective media

CRISPRi treated conditions get Dox/TMP/Media/Puro mix

CRISPRa treated conditions get Dox/Media/Puro mix

CRISPRi/a untreated conditions get normal media/puro mix

8. Plate cells (2 mL per well) and incubate

Day 2

9. Transfer cells to 1.5 mL Eppendorf tubes and spin down at 400 xg for 5 minutes.

10. Make working dilutions of Dox and TMP

Add 1 uL of Dox stock (50 mg/mL) to 999 uL media/puro

Ad 10 uL of TMP stock (100 mM) to 990 uL media/puro

11. Make condition-specific media for treated conditions

Add 15 mL of media/puro with 92.4 uL of WD Dox and 75 uL WD TMP

Add 15 ml of media/puro with 92.4 uL of WD Dox

12. Remove supernatant and resuspend with 1 mL of condition-specific media

13. Take half (0.5 mL) and plate in new well plate

14. Add 1 more mL of condition-specific media

Day 5

15. Transfer cells to 15 mL falcon tubes and spin down at 400 xg for 5 minutes.

16. Wash: Resuspend in 200 uL PBS, spin down at 300 xg for 5 min at 22C and remove supernatant

17. Resuspend cells in 350 uL RLT buffer

18. Transfer to Eppendorf tubes and freeze -80C (add to dry ice first)

--> qPCR readout